




Product Datasheet

IRF3 antibody (orb349217)

Description	Mouse monoclonal antibody to IRF3.
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Target	IRF3
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	47219
Uniprot ID	Q14653
Dilution Range	IF/ICC: 1:25, FACS: 1:25, IHC-P: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human kidney...



Flow cytometric analysis of Jurkat cells...