




Product Datasheet

SLC22A2 antibody (orb349229)

Description	Rabbit polyclonal antibody to SLC22A2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	SLC22A2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62581
Uniprot ID	O15244
Dilution Range	WB: 1:2000, IHC-P: 1:25, IHC-P: 1:25, FC: 1:25
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human stomach cancer tissue.



Immunohistochemical staining of human kidney cancer tissue.