

Product Datasheet

TUSC2 antibody (orb349286)

Description	Mouse monoclonal antibody to TUSC2.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	TUSC2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	12074
Uniprot ID	075896
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

