



Product Datasheet

BMP4 antibody (orb349310)

Description	Rabbit polyclonal antibody to BMP4.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	BMP4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46555
Uniprot ID	P12644
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



WB analysis
of A549
whole cell
lysate (L...