



Product Datasheet

CYP2E1 antibody (orb349312)

Description	Rabbit polyclonal antibody to CYP2E1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	CYP2E1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56849
Uniprot ID	P05181
Dilution Range	WB: 1:4000
Expiration Date	12 months from date of receipt.



WB analysis
of human
liver lysate
(Lane1...