

Product Datasheet

FBXO33 antibody (orb349313)

Description	Rabbit polyclonal antibody to FBXO33.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	FBXO33
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62630
Uniprot ID	Q7Z6M2
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

