




Product Datasheet

APEX1 antibody (orb349319)

Description	Mouse monoclonal antibody to APEX1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	APEX1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	35555
Uniprot ID	P27695
Dilution Range	IF/ICC: 1:25, IHC-P: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human st...



Immunofluorescence analysis of Human Cer...