



Product Datasheet

CDK5 antibody (orb349322)


Description	Mouse monoclonal antibody to CDK5.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Target	CDK5
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	33304
Uniprot ID	Q00535
Dilution Range	IF/ICC: 1:25, FACS: 1:25, IHC-P: 1:25, WB: 1:2000
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of K562 cells
U...



Immunofluorescence analysis of A549 cell...



Immunohistochemical staining of human sk...