

Product Datasheet

CRIPAK antibody (orb349346)

Description	Rabbit polyclonal antibody to CRIPAK.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	This CRIPAK antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 134-165 amino acids from the Central region of human CRIPAK.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q8N1N5
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

