



Product Datasheet

ERVK-7 antibody (orb349387)

Description	Rabbit polyclonal antibody to ERVK-7.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	ERVK-7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	66649
Uniprot ID	P61567
Dilution Range	IHC-P: 1:25, WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pl...