

## Product Datasheet

### RAB27B antibody (orb349467)

|                        |                                                                                                                                                                     |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Mouse monoclonal antibody to RAB27B.                                                                                                                                |
| <b>Species/Host</b>    | Mouse                                                                                                                                                               |
| <b>Reactivity</b>      | Human, Mouse                                                                                                                                                        |
| <b>Conjugation</b>     | Unconjugated                                                                                                                                                        |
| <b>Target</b>          | RAB27B                                                                                                                                                              |
| <b>Form/Appearance</b> | Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS. |
| <b>Note</b>            | For research use only                                                                                                                                               |
| <b>Isotype</b>         | IgG1,k                                                                                                                                                              |
| <b>Clonality</b>       | Monoclonal                                                                                                                                                          |
| <b>MW</b>              | 24608                                                                                                                                                               |
| <b>Uniprot ID</b>      | <b>O00194</b>                                                                                                                                                       |
| <b>Dilution Range</b>  | WB: 1:2000                                                                                                                                                          |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                                                     |

