



Product Datasheet

KCTD11 antibody (orb349470)

Description	Rabbit polyclonal antibody to KCTD11.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthesized Peptide
Target	KCTD11
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25887
Uniprot ID	Q693B1
Dilution Range	WB: 1:1000-1:2000, IHC-P: 1:25
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ki...