



Product Datasheet

GNG2 antibody (orb349478)

Description	Rabbit polyclonal antibody to GNG2.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthesized Peptide
Target	GNG2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	7850
Uniprot ID	P59768
Dilution Range	IHC-P: 1:25, WB: 1:8000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human br...