



Product Datasheet

NMNAT1 antibody (orb349497)

Description	Rabbit polyclonal antibody to NMNAT1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	NMNAT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31932
Uniprot ID	Q9HAN9
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis
of human
brain lysate
using ...