




Product Datasheet

RAB5B antibody (orb349522)

Description	Mouse monoclonal antibody to RAB5B.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	RAB5B
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	23707
Uniprot ID	P61020
Dilution Range	WB: 1:2000, WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis
of HeLa
whole cell
lysate (L...



WB analysis
of A431
whole cell
lysate us...