




Product Datasheet

Cleaved-Caspase 3 antibody (orb34956)

Description	Rabbit polyclonal antibody to Cleaved-Caspase 3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, IHC-P
Immunogen	Synthetic Peptide
Target	CASP3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Concentration	Lot specific, refer to datasheet supplied with delivered product.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31608
Uniprot ID	P42574
NCBI	NP_116786.1 , NP_004337.2
Dilution Range	DOT: 1:500, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Dot blot analysis of nitrocellulose memb...



Western blot analysis of Cleaved-CASP3 (...)