

Product Datasheet

NRCAM antibody (orb349573)

Description	Rabbit polyclonal antibody to NRCAM.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	This NRCAM antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide between 30-62 amino acids of mouse NRCAM.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	Q92823, Q810U4
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

