



Product Datasheet

XRN1 antibody (orb349582)

Description	Mouse monoclonal antibody to XRN1.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	XRN1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	194107
Uniprot ID	Q8IZH2
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis
of K562
whole cell
lysate (L...