

Product Datasheet

SRC antibody (orb349587)

Description	Mouse monoclonal antibody to SRC.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	SRC
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2b,k
Clonality	Monoclonal
MW	59835
Uniprot ID	P12931
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

