

## Product Datasheet

### KCTD6 antibody (FITC) (orb34960)

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>     | KCTD6 antibody (FITC)                                                       |
| <b>Species/Host</b>    | Rabbit                                                                      |
| <b>Reactivity</b>      | Human                                                                       |
| <b>Conjugation</b>     | FITC                                                                        |
| <b>Immunogen</b>       | Recombinant Human BTB/POZ domain-containing protein KCTD6 protein (1-237AA) |
| <b>Target</b>          | KCTD6                                                                       |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol              |
| <b>Form/Appearance</b> | Liquid                                                                      |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.               |
| <b>Note</b>            | For research use only                                                       |
| <b>Isotype</b>         | IgG                                                                         |
| <b>Clonality</b>       | Polyclonal                                                                  |
| <b>Uniprot ID</b>      | <a href="#">Q8NC69</a>                                                      |
| <b>Dilution Range</b>  | ELISA                                                                       |
| <b>Expiration Date</b> | 12 months from date of receipt.                                             |