



Product Datasheet

ADAMTS17 antibody (orb349663)

Description	Rabbit polyclonal antibody to ADAMTS17.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	ADAMTS17
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	121127
Uniprot ID	Q8TE56
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



WB analysis
of human
brain lysate
using ...