

Product Datasheet

SLC29A2 antibody (orb349707)

Description	Rabbit polyclonal antibody to SLC29A2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthesized Peptide
Target	SLC29A2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	50113
Uniprot ID	Q14542
Dilution Range	WB: 1:2000, IHC-P: 1:25, IHC-P: 1:25, FC: 1:25
Expiration Date	12 months from date of receipt.

