

Product Datasheet

APG8b (phospho-Thr93/Tyr99) antibody (orb34971)

Description	Rabbit polyclonal antibody to APG8b (phospho-Thr93/Tyr99)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	MAP1LC3B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14688
Uniprot ID	Q9GZQ8
NCBI	NP_073729.1
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.

