

Product Datasheet

ATP5G1 antibody (orb349713)

Description	Rabbit polyclonal antibody to ATP5G1.
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	ATP5MC1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14277
Uniprot ID	P05496
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

