

Product Datasheet

PGK1 antibody (orb349721)

Description	Mouse monoclonal antibody to PGK1.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant Protein
Target	PGK1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2a,k
Clonality	Monoclonal
MW	44615
Uniprot ID	P00558
Dilution Range	WB: 1:8000
Expiration Date	12 months from date of receipt.

