

Product Datasheet

PFN1 antibody (orb349725)

Description	Rabbit polyclonal antibody to PFN1.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthesized Peptide
Target	PFN1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	15054
Uniprot ID	P07737
Dilution Range	WB: 1:2000, FC: 1:25, FC: 1:25
Expiration Date	12 months from date of receipt.

