

Product Datasheet

RAB20 antibody (orb349736)

Description	Mouse monoclonal antibody to RAB20.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthesized Peptide
Target	RAB20
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	26277
Uniprot ID	Q9NX57
Dilution Range	WB: 1:2000-1:4000, FC: 1:25
Expiration Date	12 months from date of receipt.

