



Product Datasheet

Cox4i1 antibody (orb349737)

Description	Rabbit polyclonal antibody to Cox4i1.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthesized Peptide
Target	Cox4i1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19530
Uniprot ID	P19783
Dilution Range	WB: 1:2000, FC: 1:25
Expiration Date	12 months from date of receipt.



WB analysis of NIH/3T3 whole cell lysate...