

Product Datasheet

CSRNP1 antibody (orb349760)

Description	Rabbit polyclonal antibody to CSRNP1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	CSRNP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	63508
Uniprot ID	Q96S65
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

