

Product Datasheet

COX IV antibody (orb349762)

Description	Rabbit polyclonal antibody to COX IV
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	COX4I2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20010
Uniprot ID	Q96KJ9
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

