

Product Datasheet

S1 tag antibody (orb349764)

Description	Rabbit polyclonal antibody to S1 tag.
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	This S1 tag antibody is generated from a rabbit immunized with a KLH conjugated synthetic peptide of S1 tag.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	NA
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

