

Product Datasheet

BCDIN3D antibody (orb349775)

Description	Rabbit polyclonal antibody to BCDIN3D.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	BCDIN3D
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33200
Uniprot ID	Q7Z5W3
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.

