

Product Datasheet

DBN1 antibody (orb349792)

Description	Rabbit polyclonal antibody to DBN1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	DBN1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	71429
Uniprot ID	Q16643
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.

