



Product Datasheet

Sox2 (phospho-Ser246) antibody (orb34982)

Description	Rabbit polyclonal to Sox2 (phospho-Ser246)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, IF
Immunogen	Synthetic Peptide
Target	SOX2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	34310
Uniprot ID	P48431
NCBI	NP_003097.1
Dilution Range	IF/ICC: 1:200, DOT: 1:500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SY5Y cell...