



Product Datasheet

GATA6 (phospho-Tyr271) antibody (orb34996)

Description	Rabbit polyclonal to GATA6 (phospho-Tyr271)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, IF
Immunogen	Synthetic Peptide
Target	GATA6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60033
Uniprot ID	Q92908
NCBI	NP_005248.2
Dilution Range	IF/ICC: 1:100, DOT: 1:500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cells