



Product Datasheet

PTCD2 antibody (orb35013)

Description	PTCD2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Pentatricopeptide repeat-containing protein 2, mitochondrial protein (153-388AA)
Target	PTCD2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WV60
Dilution Range	WB:1:2000-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse heart & k...



Immunohistochemical staining of human pr...