



Product Datasheet

ACAA1 antibody (orb350245)

Description	ACAA1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human 3-ketoacyl-CoA thiolase, peroxisomal protein (27-331AA)
Target	Disintegrin EC3B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09110
Dilution Range	WB:1:500-5000, IHC-P:1:200-500, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human th...



Immunohistochemical staining of human pr...