

Product Datasheet

ADAMTS4 antibody (FITC) (orb350255)

Description	ADAMTS4 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human A disintegrin and metalloproteinase with thrombospondin motifs 4 protein (213-319AA)
Target	Prothoracicotropic hormone
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075173
Dilution Range	IHC-P:1:500-1000, IF:1:50-200
Expiration Date	12 months from date of receipt.