

Product Datasheet

Annexin A5 antibody (orb350277)

Description	Annexin A5 antibody
Species/Host	Rabbit
Reactivity	Reptile
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Cynops pyrrhogaster Annexin A5 protein (1-323AA)
Target	Annexin A5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P70075
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.