



Product Datasheet

BMPR1A antibody (orb350325)


Description	BMPR1A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Bone morphogenetic protein receptor type-1A protein (177-532AA)
Target	BMPR1A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P36894
Dilution Range	WB:1:500-5000, IHC-P:1:20-1000, IF:1:200-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human kidney tissue.



Immunofluorescence analysis of HepG2 cells.



Immunohistochemical staining of human placenta tissue.