



Product Datasheet

CDK5 antibody (orb350357)


Description	CDK5 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Cyclin-dependent-like kinase 5 protein (1-292AA)
Target	CDK5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q00535
Dilution Range	WB:1:200-1000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse brain tis...



Western blot analysis of mouse heart & b...



Immunohistochemical staining of human ov...