

Product Datasheet

CNTFR antibody (FITC) (orb350371)

Description	CNTFR antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Ciliary neurotrophic factor receptor subunit alpha protein (106-342AA)
Target	CNTFR
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P26992
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.