



Product Datasheet

CXCR4 antibody (orb350401)


Description	CXCR4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human C-X-C chemokine receptor type 4 protein (1-356AA)
Target	CXCR4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P61073
Dilution Range	WB:1:500-2000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela cells using CXCR4 antibody (orb350401).



Western blot analysis of mouse liver & brain tissue using CXCR4 antibody (orb350401).



Immunofluorescence analysis of Hela cells using CXCR4 antibody (orb350401).