



Product Datasheet

DDX3X antibody (orb350425)


Description	DDX3X antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human ATP-dependent RNA helicase DDX3X protein (2-662AA)
Target	DDX3X
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	000571
Dilution Range	WB:1:500-5000, IHC-P:1:100-500, IF:1:50-500, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH3T3(Lane 1),...



Western blot analysis of K562(Lane 1), J...



Immunofluorescence analysis of HepG2 cel...