



Product Datasheet

CQ068 antibody (orb35057)

Description	Rabbit polyclonal antibody to CQ068. The protein is a member of the CTC1 family, and it is subunits o...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	CTC1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	134609
Uniprot ID	Q2NKJ3
NCBI	NP_079375.3
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
CQ068
antibo...