



Product Datasheet

MAPK9 antibody (orb350582)


Description	MAPK9 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP
Immunogen	Recombinant Human Mitogen-activated protein kinase 9 protein (1-424AA)
Target	MAPK9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P45984
Dilution Range	IF:1:200-500, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...



Immunoprecipitation analysis of HeLa who...



Immunofluorescence analysis of HeLa cell...