

Product Datasheet

S1PR1 antibody (FITC) (orb350687)

Description	S1PR1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Sphingosine 1-phosphate receptor 1 protein (311-382AA)
Target	S1PR1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P21453
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.