

Product Datasheet

Alpha-2 antiplasmin antibody (FITC) (orb350698)

Description	Alpha-2 antiplasmin antibody (FITC)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	FITC
Immunogen	Recombinant Rat Alpha-2 antiplasmin protein (20-418AA)
Target	Serpinf1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q80ZA3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.