

Product Datasheet

TRIM43 antibody (HRP) (orb35071)

Description	TRIM43 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Tripartite motif-containing protein 43 protein (43-240AA)
Target	TRIM43
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96BQ3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.