



Product Datasheet

SOD2 antibody (orb350712)


Description	SOD2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Superoxide dismutase [Mn], mitochondrial protein (64-175AA)
Target	SOD2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04179
Dilution Range	WB:1:500-5000, IHC-P:1:100-600, IF:1:200-500, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Hela whole cel...



Western blot aanalysis of HepG2 whole ce...